

# SORFA® SAMPLING SYSTEM

## NASAL SWAB SAMPLE COLLECTION KIT



**SPARTAN**<sup>™</sup>  
M E D I C A L

## RELIABLE SOLUTION FOR SAMPLE COLLECTION OF COVID-19, INFLUENZA AND OTHER VIRUSES

### Stability under a Wide Temperature Range

The systems can be flexibly transported and utilized under a wide temperature range from -70°C to 25°C; Can be stored between 15-20°C

### Usability for Safe & Convenient Viral Testing

State-of-the-art swabs, outstanding medium formulation and leak-proof transport tubes deliver great reliability

### Reasonable Structure

- Superior leak-proof tubes with conical bottom ideal sample centrifugation
- Nasopharyngeal swabs are easily fractured, and armed with efficient flocked tips for sampling
- Two beads in each tube designed for rapid specimen release and dispersion
- Ergonomic structure of swabs refers to facial structural data of different people, delivers safe and efficient sampling process



### KIT COMPONENTS:


- One Flocked Nasal Swab
- Transfer Media (2 or 3mL Liquid)
- Test Tube (5 or 10mL)

### VTM-N Sampling System

- ✓ VTM-N medium for transportation and prevention of inactive virus, inactive microbial viability enables safe handling, transport and processing.

### VTM Sampling System

- ✓ The VTM medium is preferably applied for transportation and preservation of active virus
- ✓ Prefilled with PH indicator for visible judgement of pollution



PART NO.	PRODUCT	DESCRIPTION
751410	SORFA Sampling System - 2ml Viral Transport Media (VTM)	Flocking Nasopharyngeal Swab, 5mL Tube with 2mL Viral Transport Media (VTM) - 50 Swab Test in Box
751420	SORFA Sampling System - 3ml Viral Transport Media (VTM)	Flocking Nasopharyngeal Swab, 10ml tube with 3ml Viral Transport Media (VTM) - 50 Swab Test in Box



### LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

[CS@SPARTANMEDSPINE.COM](mailto:CS@SPARTANMEDSPINE.COM)



### SPARTAN MEDICAL INC.

1445 Research Blvd, Suite 520,  
Rockville, MD 20850

[spartanmedspine.com](http://spartanmedspine.com)

DUNS: 807438341 | CAGE:65VT5

